



Safety Check & Tracking Sheet

Always inspect used equipment before taking it home. Ask for instructions and maintenance tips.

Quick Equipment Safety Check

- All parts present (no missing bolts/pins).
- No cracks, sharp edges, or bent frame parts.
- Brakes work and lock securely.
- Wheels/tires roll smoothly (no wobble); casters spin freely.
- Seat/backrest stable; armrests/footrests secure (wheelchair).
- Height adjustments lock firmly (walker/cane).
- No frayed straps or damaged slings (patient lift).
- Electrical cords/battery connections intact (power devices).
- I received basic instructions and know who to call for help/repairs.

Program Tracker

Program / Supplier	Contact	Item requested	Status (date & next step)